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Implantologie

Implantology

125

Implantologie
Implantology

3. Implantologie · Implantology

Knochenmühle · Bone mill

1720

8,0 x 12,9 cm
Knochenmühle
Bone mill

126

i

Knochenmühle · Bone mill

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Bone Mill for the preparation of autologous bone material. In order to correct extensive bone defects, an intraoral extraction of bone material from another site of the jaw is necessary in many cases. Larger bone pieces or bone blocks are extracted with the help of e.g. a trephine drill.

The use of the bone mill enables an effective milling of the bone material in order to produce a volume of bone grafts as large as possible.

The “granular” consistency of the produced bone graft ensures an optimal adaptation to the bone base.

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Bone mill forceps for the preparation of autologous bone material.

- Exchangeable titanium bone mill inserts
- Easy to clean
- No loss of bone material
- Conform to the RKI-guidelines
- Specially designed new forceps, which allows parallel closing of the milling parts
- Maximum milling result
- Stop screw avoids direct contact of the milling parts and protects from wear

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Insert the bone piece into the crusher and use the pestle in order to crush the bone. If necessary, you may use a mallet additionally. The SCHWERT bone crusher is a very time- and cost efficient alternative for the bone mill 1720.

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
The bone scraper is intended for gaining of autologous bone material. This manual working assures an optimal cell viability and reduced absorption. The bone is collected in the underlying chamber while scraping and mixed up with blood simultaneously. The bone material may be used directly out of the chamber. The blade is suited for up to 30 applications.

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With the bone scraper by Buser is an instrument provided which has working tips with the measures of 3 mm and 5 mm. This instrument makes it very easy to collect autologous bone directly in the area of the treatment. Specially hardened working tips ensure sharpness for a long time.

Beside that the instruments suits great for contouring the alveolar ridge.

129

1418-3

3 mm / 5 mm

2616-25

Medizinbecher, 25 ccm,
rostfrei
Medicine cup, 25 ccm,
stainless steel

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Bone Trap / Suction tube incl. filter. Already during the preparation of the implant site, the bone trap may be used to gain autologous bone material. It is completely dismountable for cleaning and maintenance. The stainless steel filter is designed for multiple use.

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
Osteotomes are used for increasing the initial stability of the implant. The stopping screw increases the safety during the sinus lift significantly.

Artikel-Nr. Item-No.	Form der Instrumentenspitze Shape of the instrument tip	Ø außen Ø exterior	Graduierung graduation
1703-127	konkav, konisch, bajonett concave, conical, bayonet	2,7 mm	7 / 10 / 13 / 15 / 18 / 20
1703-132	konkav, konisch, bajonett concave, conical, bayonet	3,2 mm	7 / 10 / 13 / 15 / 18 / 20
1703-137	konkav, konisch, bajonett concave, conical, bayonet	3,7 mm	7 / 10 / 13 / 15 / 18 / 20
1703-142	konkav, konisch, bajonett concave, conical, bayonet	4,2 mm	7 / 10 / 13 / 15 / 18 / 20
1703-150	konkav, konisch, bajonett concave, conical, bayonet	5,0 mm	7 / 10 / 13 / 15 / 18 / 20

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

Application: After the pilot drill is carried out, the first instrument is driven in by hand or with gently mallet blows (1637 / Page 82). The instrument tip is adjusted parallel to the alveolar crest. After driving it in up to the shank end, careful turning movement, are performed for bone spreading. The turning has to be performed maximally until it is possible to insert the next larger instrument. The steps are repeated, up to the required diameter for the implant is reached. The gap can be rounded out according to the implant by turning movements of the instrument's handle.

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
An optimum sensitivity during the laterally lift procedure.

The perforation of the Schneiderian membrane is a complication that often occurs during a maxillary sinus floor augmentation procedure. The new SCHWERT Sinus-Lift Instruments increase the sensitivity by blunt and flexible working tips with "Swing-Effect". Because of their special shapes and the characteristic of the stainless steel between the handle and the active part, light pressure is cushioned and gives the operator a direct and sensitive response. The risk to perforate the Schneiderian membrane is reduced to a minimum. Already during their development the focus was a practicable handling. Because of varying angulations and shapes of the active parts they are easy to navigate also in difficult areas. The line of currettes consists of five different figures for every situation. Spoon and condenser complement the set to place the bone graft.

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


SCHWERT cares!
From each instrument sold 1.50 Euro goes to a fund for disadvantaged children.

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1. The tip of the instrument serves to determine the depth of the osteotomy. Scaling at the tip of 8 / 10 / 11,5 / 13 / 15 / 17 mm in order to place the implant into the bone. The diameter of the tip is 2,2 mm for the positioning of the first drilling perforation.

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2. The length of the sliding caliper enables the measuring and length check of the cutter for the first drill.

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3. The rectangular part at the distal end of the instrument has a horizontal as well as a vertical measurement. (**1734-1** 7 mm • **1734-2** 9 mm)

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4. The horizontal measurement up to 7 mm / 9 mm enables the positioning of an implant of 3,4 mm - 6 mm diameter – concerning the distance between the implants.

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5. The vertical measurement up to 8 mm enables the positioning of the implant of a diameter of 4-7 (4-8 mm) mm or 5-8 mm (5-9 mm).this has to be repeated on the opposite side of the tooth.

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6. If there is a change of positioning the implant it is possible to check 2 mm space to prevent a perforation of the interproximale bone.

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The bone block is adapted with the curved forceps. The block is safely fixed, from the first drill to the final fixation with a screw.

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Parallel drilling guide, for easy parallel implant insertion. (Ø 2.1 for pilot drill, 1.85 mm size of guiding pin) insert the parallel drilling guide into the first pilot drill with the guiding pin. With the fixable joint a precise adjustment for the following parallel pilot drilling is possible.

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
Inserting Implants gets much easier and effective due the caliper acc. to Prof. Dr. med. dent. Gisbert Krekeler. Working is safer, faster and more precise by the combination several functions in one instrument.

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Functionality of the Krekeler Caliper

For a maximum stability always the widest possible implant diameter should be chosen. Due to the exactly 1,5 mm measuring rods, you can read the maximum diameter directly. Due to scientific research in order to avoid bones resorbition after the implantological treatment, there should be a wall of at least 1,5 mm of bone material to each side.

1. During the measurement of the interdental gap, it is possible to read directly the maximum implant diameter from the IN mark. A correct measurement is performed by applying the rods at the cervix dentis of the neighboring teeth. From the OUT mark you may read the the overall measure.

2. The sharpened rods also support the operator in placing the first mark in the bone, to ensure a safe distance to the neighbouring implant / tooth. If there is a neighbour implant it is possible to use one rod as spacer for the minimum distance of the next drill.

3. At the distal ending of the Krekeler Caliper you also find a drilling guide and a 5 and 7 mm graduation to measure the patients alveolar ridge from crestal. This supports an easy and direct evaluation if there is enough bone material left.

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
Acc. to Dr. J.C.M. Rosa. These selection of the right chisel to remove the graft from maxillary tuberosity is performed according to the shape of the existing defect in the receiving area and required bone volume. The most suitable chisel shapes are the straight and the gouge ones of 6, 8 or 10 mm width; as a general rule, the chosen chisel should be 2 mm wider than the graft width to be removed. These chisels have a millimetered active part with the aim to ease the graft measurement at the moment of removal, allowing the obtention of a bone specimen with suitable length and thickness.

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
Instrument kit according Dr. Rosa, prepared for the surgical protocol for the IDR technique.

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