


ABACUS® Diamantinstrumente • ABACUS® Diamond instruments • Instruments diamantés ABACUS®

Z847KR (546)

Konus, Kante rund
KR taper, modified shoulder
Cône, bout arrondi



Größe/Size/Taille ø 1/10 mm	016	018	
Länge/Length/Longueur mm	8,0	8,0	
Winkel/Angle/Angle α	1,8°	1,6°	

5	REF			
	Z847KR-	016	018	C-FG

Z848 (173)

Konus, flach
Flat end taper
Cône, bout plat



Größe/Size/Taille ø 1/10 mm	018		
Länge/Length/Longueur mm	10,0		
Winkel/Angle/Angle α	2,0°		

5	REF			
	Z848-	018	C-FG	

Z850 (199)

Konus, rund
Round end taper
Cône, bout arrondi



Größe/Size/Taille ø 1/10 mm	012	014	016	
Länge/Length/Longueur mm	10,0	10,0	10,0	
Winkel/Angle/Angle α	1,9	1,7	2,0	

5	REF				
	Z850-	012	014	016	C-FG
	Z850-	012	014		M-FG
	Z850-	012			F-FG

Z855 (197)

Konus, rund
Round end taper
Cône, bout arrondi



Größe/Size/Taille ø 1/10 mm	025		
Länge/Length/Longueur mm	7,0		
Winkel/Angle/Angle α	4,5°		

5	REF			
	Z855-	025	C-FG	

Z856 (198)

Konus, rund
Round end taper
Cône, bout arrondi



Größe/Size/Taille ø 1/10 mm	014	016	018	021	
Länge/Length/Longueur mm	8,0	8,0	8,0	8,0	
Winkel/Angle/Angle α	1,5°	1,4°	1,4°	2,8°	

5	REF					
	Z856-	014	016	018	021	C-FG

Z856L (201)

Konus, rund, lang
Long round end taper
Cône rond



Größe/Size/Taille ø 1/10 mm	018		
Länge/Length/Longueur mm	9,0		
Winkel/Angle/Angle α	1,4°		

5	REF			
	Z856L-	018	C-FG	

Z858 (165)

Konus, spitz
Needle
Cône pointu



Größe/Size/Taille ø 1/10 mm	014		
Länge/Length/Longueur mm	8,0		
Winkel/Angle/Angle α	3,0°		

5	REF			
	Z858-	014	C-FG	

Z862 (249)

Flamme
Flame
Flamme



Größe/Size/Taille ø 1/10 mm	014		
Länge/Length/Longueur mm	8,0		

5	REF			
	Z862-	014	C-FG	

Z863 (250)

Flamme
Flame
Flamme



Größe/Size/Taille ø 1/10 mm	012	016	
Länge/Length/Longueur mm	10,0	10,0	

5	REF			
	Z863-	012	016	C-FG
	Z863-	012	016	M-FG
	Z863-	012	016	F-FG

Z877K (297)

Torpedo, konisch
Torpedo, conical
Torpille conique



Größe/Size/Taille ø 1/10 mm	012	016	
Länge/Length/Longueur mm	6,0	6,0	
Winkel/Angle/Angle α	2,1°	3,0°	

5	REF			
	Z877K-	012	016	C-FG
	Z877K-		016	F-FG

Bitte Drehzahlbereiche beachten (Drehzahltable - Seite 9.2) / Please observe recommended speeds (speed chart - page 9.2) / Respecter les vitesses de rotation recommandées (voir à la page 9.2)





	... C	Grüner Ring/Green ring/Bague verte	Grob/Coarse/Grain gros	175 µm	ISO No. 534
	... M	Blauer Ring/Blue ring/Bague bleue	Mittel/Medium/Grain moyen	125 µm	ISO No. 524
	... F	Roter Ring/Red ring/Bague rouge	Fein/Fine/Grain fin	45 µm	ISO No. 514

**Z878
(289)**

Torpedo
Beveled cylinder
Torpille




Größe/Size/Taille ø 1/10 mm	010	012	014
Länge/Length/Longueur mm	8,0	8,0	8,0





 5	REF			
	Z878-	012	014	C-FG
	Z878-	010	012	M-FG
	Z878-	010	012	F-FG

**Z878K
(298)**

Torpedo, konisch
Torpedo, conical
Torpille conique



Größe/Size/Taille ø 1/10 mm	014	016	018	020	022
Länge/Length/Longueur mm	8,0	8,0	8,0	8,0	8,0
Winkel/Angle/Angle α	2,0°	2,1°	2,5°	2,5°	2,5°





 5	REF					
	Z878K-	014	016	018	020	C-FG
	Z878K-	014	016			M-FG
	Z878K-		016		022	F-FG

**Z879
(290)**

Torpedo
Beveled cylinder
Torpille



Größe/Size/Taille ø 1/10 mm	012	014
Länge/Length/Longueur mm	10,0	10,0





 5	REF			
	Z879-	012	014	C-FG
	Z879-	012	014	M-FG
	Z879-	012	014	F-FG

**Z879K
(299)**

Torpedo, konisch
Torpedo, conical
Torpille conique



Größe/Size/Taille ø 1/10 mm	012	014	016	018
Länge/Length/Longueur mm	10,0	10,0	10,0	10,0
Winkel/Angle/Angle α	1,6°	1,6°	1,6°	1,8°




 5	REF					
	Z879K-	012	014	016	018	C-FG
	Z879K-	012	014	016		M-FG
	Z879K-	012	014	016		F-FG

**Z880
(140)**

Zylinder, rund
Round end cylinder
Cylindre, bout arrondi



Größe/Size/Taille ø 1/10 mm	012
Länge/Length/Longueur mm	6,0





 5	REF		
	Z880-	012	C-FG
	Z880-	012	M-FG

**Z881
(141)**

Zylinder, rund
Round end cylinder
Cylindre, bout arrondi



Größe/Size/Taille ø 1/10 mm	014	016
Länge/Length/Longueur mm	8,0	8,0

 5	REF			
	Z881-	014	016	C-FG
	Z881-	014		M-FG
	Z881-	014	016	F-FG

**Z909
(068)**

Reifen
Donut
Roue






Größe/Size/Taille ø 1/10 mm	040
Länge/Length/Longueur mm	1,6

 5	REF		
	Z909-	040	C-FG



Bitte Drehzahlbereiche beachten (Drehzahltable - Seite 9.2) / Please observe recommended speeds (speed chart - page 9.2) / Respecter les vitesses de rotation recommandées (voir à la page 9.2)

	... C	Grüner Ring/Green ring/Bague verte	Grob/Coarse/Grain gros	175 µm	ISO No. 534
	... M	Blauer Ring/Blue ring/Bague bleue	Mittel/Medium/Grain moyen	125 µm	ISO No. 524
	... F	Roter Ring/Red ring/Bague rouge	Fein/Fine/Grain fin	45 µm	ISO No. 514