

PROTHÈSE

PORTE-EMPREINTE

PROTHÈSE

PORTE-EMPREINTE


N°2 BORDS RÉTENTIFS

SUPÉRIEUR




IDEAL POUR ALGINATE

INFÉRIEUR



Taille mm	REF	Taille standard
61x47	3041.U1	1
64x51	3041.U2	2
66x53	3041.U3	3
69x57	3041.U4	4
71x60	3041.U5	5
73x62	3041.U6	6




Taille mm	REF	Taille standard
55x47	3041.L1	1
60x49	3041.L2	2
72x54	3041.L3	3
75x57	3041.L4	4
76x59	3041.L5	5
78x61	3041.L6	6

N°1 BORDS MINCES DROITS


SUPÉRIEUR



INFÉRIEUR



Taille mm	REF	Taille standard
61x47	3040.U1	1
64x51	3040.U2	2
66x53	3040.U3	3
69x57	3040.U4	4
71x60	3040.U5	5
73x62	3040.U6	6




Taille mm	REF	Taille standard
55x47	3040.L1	1
60x49	3040.L2	2
72x54	3040.L3	3
75x57	3040.L4	4
76x59	3040.L5	5
78x61	3040.L6	6

N°3 PLEIN, BORDS RÉTENTIFS


SUPÉRIEUR



INFÉRIEUR



Taille mm	REF	Taille standard
61x47	3050.U1	1
64x51	3050.U2	2
66x53	3050.U3	3
69x57	3050.U4	4
71x60	3050.U5	5
73x62	3050.U6	6



Taille mm	REF	Taille standard
55x47	3050.L1	1
60x49	3050.L2	2
72x54	3050.L3	3
75x57	3050.L4	4
76x59	3050.L5	5
78x61	3050.L6	6

DEMI PORTE-EMPREINTE

DROIT



GAUCHE



	Taille mm	REF
GAUCHE	30x49	3040.L2S
	37x57	3040.L4S
	39x61	3040.L6S
DROIT	30x49	3040.L2D
	37x57	3040.L4D
	39x61	3040.L6D

Personnalisez vos porte-empreintes gratuitement à partir de 20 pièces.